




YOU

**& BE<sup>®</sup>**

A COMPLETE SOLUTION DESIGNED FOR YOUR PATIENTS' UROSTOMY NEEDS

PRODUCT GUIDE SUPPLEMENT



**CONFIDENCE BE<sup>®</sup>**

CONFIDENCE BE <sup>®</sup> UROSTOMY							
LARGE BAGS							
SIZE	STARTER HOLE	PRE-CUT	Transparent Stone	Trans/Overlap Black	Trans/Overlap White	Trans/Overlap Stone	PACK SIZE
13mm	•	–	BULT13ST	BUL13B	BUL13W	BUL13ST	10
STANDARD BAGS							
SIZE	STARTER HOLE	PRE-CUT	Transparent Stone	Trans/Overlap Black	Trans/Overlap White	Trans/Overlap Stone	PACK SIZE
13mm	•	–	BUT13ST	BU13B	BU13W	BU13ST	10
21mm	–	•	–	BU21B	BU21W	BU21ST	10
25mm	–	•	–	BU25B	BU25W	BU25ST	10
28mm	–	•	–	BU28B	BU28W	BU28ST	10
32mm	–	•	–	BU32B	BU32W	BU32ST	10
35mm	–	•	–	BU35B	BU35W	BU35ST	10
CONFIDENCE BE <sup>®</sup> UROSTOMY SOFT CONVEX							
LARGE BAGS							
SIZE	STARTER HOLE	PRE-CUT	Transparent Stone	Trans/Overlap Black	Trans/Overlap White	Trans/Overlap Stone	PACK SIZE
13–25mm	•	–	XBULT1325ST	XBUL1325B	XBUL1325W	XBUL1325ST	10
13–38mm	•	–	XBULT1338ST	XBUL1338B	XBUL1338W	XBUL1338ST	10
13–52mm	•	–	XBULT1352ST	XBUL1352B	XBUL1352W	XBUL1352ST	10
STANDARD BAGS							
SIZE	STARTER HOLE	PRE-CUT	Transparent Stone	Trans/Overlap Black	Trans/Overlap White	Trans/Overlap Stone	PACK SIZE
13–25mm	•	–	XBUT1325ST	XBU1325B	XBU1325W	XBU1325ST	10
13–38mm	•	–	XBUT1338ST	XBU1338B	XBU1338W	XBU1338ST	10
13–52mm	•	–	XBUT1352ST	XBU1352B	XBU1352W	XBU1352ST	10
21mm	–	•	–	XBU21B	XBU21W	XBU21ST	10
25mm	–	•	–	XBU25B	XBU25W	XBU25ST	10
28mm	–	•	–	XBU28B	XBU28W	XBU28ST	10
32mm	–	•	–	XBU32B	XBU32W	XBU32ST	10
35mm	–	•	–	XBU35B	XBU35W	XBU35ST	10
CONFIDENCE BE <sup>®</sup> NIGHT DRAINAGE BAG							
SIZE	ORDER CODE						QUANTITY
2000ml	BNDB						10

For free samples or further product information, please contact your local Territory Manager or:



**FREEPHONE Nurse Direct 0800 028 8014** (England & Wales)

**FREEPHONE 0800 626388** (Scotland)

**FREEPHONE 1-800 408508** (Rep. of Ireland)



[hello@salts.co.uk](mailto:hello@salts.co.uk)



[www.salts.co.uk](http://www.salts.co.uk)

Follow us @SaltsHealthcare:    



**CE UK CA** © Registered trade marks of Salts Healthcare Ltd. © Salts Healthcare Ltd 2023. Products and trade marks of Salts Healthcare Ltd are protected by UK and foreign patents, registered designs and trade marks. VELCRO<sup>®</sup> is a registered trade mark of Velcro BVBA. Used with permission. For further details, please visit [www.salts.co.uk](http://www.salts.co.uk) **References:** 1. In comparison to Salts Confidence Natural Urostomy and Salts Confidence Convex Supersoft Urostomy bags. Data on file. Salts Healthcare 2023. 2. Confidence BE<sup>®</sup> Urostomy Phase 1 and 2 clinical surveys. Data on file. Salts Healthcare 2022. 3. Dzulkiifi, M.A., & Mustafar, M.F. (2013). The influence of colour on memory performance: a review. *The Malaysian journal of medical sciences: MJMS*, 20(2), 3–9. RM765958 03/23 v.1

# CONFIDENCE BE® UROSTOMY



Since the launch of Confidence BE® we've continued to listen to nurses and patients to design a bag to meet the needs of your urostomy patients. Available in three colours – black, white and stone – Confidence BE® Urostomy has a range of features to support health and well-being. Its soft outlet and bung is designed to be comfortable without compromising ease of use, the outlet tucks away securely for added discretion and a detachable bung lanyard can be moved away from the flow of urine when draining.

## Features:

### BE COMFORTABLE

- Soft, water-repellent fabric.
- Soft outlet and bung is comfortable without compromising ease of use.
- Urostomy hydrocolloid resists fluid and erosion.<sup>1</sup>
- Wafer is infused with aloe extracts.

### BE YOU

- The first urostomy bag to offer a choice of three colours: black, white and stone.
- Available in flat, soft convex, standard and large.

### BE DISCREET

- The sleek, low-profile design looks less like a medical device.
- Low-rustle fabric cover.<sup>2</sup>
- Outlet tucks away and is securely fastened with premium VELCRO® hook and loop fastener material.

### BE FLEXIBLE

- Unique, patented Flexifit® wafer bends and flexes with movement.

### BE IN CONTROL

- Flow control makes it easy to control the flow of urine.
- Rotating and detachable bung lanyard can be moved away from urine flow when draining.
- Quick-drying fabric.
- Blue, guided touchpoints on the outlet may make the bag easier to see and use.<sup>2,3</sup>



LARGE FLAT BAGS



STANDARD FLAT BAGS



LARGE CONVEX BAGS



STANDARD CONVEX BAGS

# CONFIDENCE BE® NIGHT DRAINAGE BAG

The Confidence BE® Urostomy bag is designed to connect to the Night Drainage Bag easily, helping to support patients day and night. Recyclable cardboard packaging for the Night Drainage Bag discreetly turns into a nightstand.

## Night Drainage Bag features:

### BE IN CONTROL

- Innovative tubing contains a filament designed to maintain flow of urine, even if kinked and compressed.
- Integrated handle to easily transport the bag for emptying.
- Visual lock indicator shows when the tap is securely in the closed position.
- Blue, guided touchpoints may make the bag easier to see and use.<sup>3</sup>
- 2000ml minimum and maximum capacity.

### BE COMFORTABLE

- Soft and flexible for more comfort, the tube (1.2m) is designed to give more freedom of movement and maintain a secure connection.

### BE YOU

- 100% PVC free for a more environmentally friendly Night Drainage Bag.

## Night Drainage Stand features:

### BE DISCREET

- Cardboard box converts into a discreet stand which obscures the contents of the Night Drainage Bag.
- Splash resistant coating makes the box more durable against liquid.

### BE IN CONTROL

- Integrated handles for easy transportation.
- Drainage hole knockout can be used to empty the bag without removing it from the stand.

### BE YOU

- Helps to reduce plastic waste from traditional night drainage stands.

